

## TG-51 Worksheet C: Kecal N<sub>D,w</sub><sup>60Co</sup> for plane-parallel chambers

There are two methods for determining Kecal N<sub>D,w</sub><sup>60Co</sup> for plane-parallel chamber. Method A uses cross-calibration against a calibrated cylindrical chamber and is preferred method. Method B uses a <sup>60</sup>Co absorbed-dose calibration factor.

### Method A: Cross-Calibration

#### 1. Site data

Institution: \_\_\_\_\_  
 Physicist: \_\_\_\_\_  
 Date: \_\_\_\_\_  
 Accel Mfr.: \_\_\_\_\_  
 Model & serial number: \_\_\_\_\_  
 Nominal e<sup>-</sup> energy/beam identifier: \_\_\_\_\_ MeV

#### 2. Dose using cylindrical chamber

Do reference dosimetry for this beam using Worksheet B.

Transfer the following information from that worksheet:

- a. Date: \_\_\_\_\_ [B:1]
- b. Beam quality  $R_{50}$ : \_\_\_\_\_ [B:4a]
- c. Reference depth,  $d_{ref}$ : \_\_\_\_\_ [B:4b]
- d. Dose/MU at  $d_{ref}$ : \_\_\_\_\_ [B:10c]
- e. Number of MU (same used here): \_\_\_\_\_ [B:3c]

Now place the point of measurement of the plane-parallel chamber at  $d_{ref}$

#### 3. Temperature/ Pressure Correction (Sec. VII. C)

a. Temperature: \_\_\_\_\_ °C

b. Pressure: \_\_\_\_\_ kPa  $= \text{mmHg} \cdot \frac{101.33}{760}$   
 $\left[ \text{Eq.(10)} = \left( \frac{273.2 + 3a}{295.2} \right) \left( \frac{101.33}{3b} \right) \right]$

#### 4. Polarity Correction (Sec. VII. A.)

$M_{raw}^+$ : \_\_\_\_\_ C or rdg

$M_{raw}^-$ : \_\_\_\_\_ C or rdg

a.  $M_{raw}$  (for polarity used clinically): \_\_\_\_\_ C or rdg

b.  $P_{pol}$ : \_\_\_\_\_  $\left[ \text{Eq.(9)} = \left| \frac{(M_{raw}^+ - M_{raw}^-)}{2M_{raw}} \right| \right]$

#### 5. $P_{ion}$ measurements (Sec. VII. D.2)

Operating Voltage =  $V_H$ : \_\_\_\_\_ V

Lower voltage  $V_L$ : \_\_\_\_\_ V

$M_{raw}^H$ : \_\_\_\_\_ C or rdg

$M_{raw}^L$ : \_\_\_\_\_ C or rdg

$P_{ion}$  ( $V_H$ ) (pulsed/swept beam, Eq.(12)): \_\_\_\_\_  $\left[ \left( 1 - \frac{V_H}{V_L} \right) \middle/ \left( \frac{M_{raw}^H}{M_{raw}^L} - \frac{V_H}{V_L} \right) \right]$

If  $P_{ion} > 1.05$ , another ion chamber should be used.

# TG-51 Worksheet C: Kecal $N_{D,w}^{60_{Co}}$ for plane-parallel chambers (cont)

## 6. Corrected ion. Ch. Rdg. M (Sec. VII)

$$M = P_{ion} P_{TP} P_{elec} (=1.0) P_{pol} M_{raw} = [5 \cdot 3c \cdot 1.0 \cdot 4b \cdot 4a]$$

Fully corrected M(Eq.(8)): \_\_\_\_\_ C or rdg

## 7. Determination of $k_{R_{50}}$ for plane-parallel chamber, beam quality $R_{50}$ (2b)

i.  $k_{R_{50}}$  from figures \_\_\_\_\_ [Fig 6 or 8]

or ii.  $k_{R_{50}}$  from analytic expression for well-guarded plane-parallel chambers

$$k_{R_{50}} = 1.2239 - 0.145(R_{50})^{0.214} \quad \text{[Eq.(20) } 2 \leq R_{50} \leq 20 \text{ cm]}$$

## 8. Cross calibration value

$$(k_{ecal} N_{D,w}^{60_{Co}})^{pp} = \frac{(D_w / MU)^{cyl} MU}{(M k_{R_{50}})^{pp}} = \left[ \frac{2d \cdot 2e}{6 \cdot 7(i \text{ or } ii)} \right] \quad \boxed{\quad \quad \quad \text{Gy/C(or Gy/rdg)}}$$

## Method B: $^{60}\text{Co}$ Calibration

### 1. Instrumentation

a. Chamber model: \_\_\_\_\_

Serial number: \_\_\_\_\_

Waterproof: yes

no

If no, is waterproofing  $\leq$  1 mm PMMA or thin latex?: yes  no

b. Electrometer model: \_\_\_\_\_

Serial number: \_\_\_\_\_

i.  $P_{elec}$ , electrom. Corr factor (Sec. VII. B): \_\_\_\_\_ C/C or C/rdg.

c. Calibration Factor  $N_{D,w}^{60_{Co}}$  (Sec. V): \_\_\_\_\_ Gy/C (or Gy/rdg)

Date of report (not to exceed 2 years): \_\_\_\_\_

### 2. Determination of $k_{ecal}$

Chamber model used to get  $k_{ecal}$ : \_\_\_\_\_

a.  $k_{ecal}$ : \_\_\_\_\_ [Table II]

### 3. $k_{ecal} N_{D,w}^{60_{Co}}$ :

Gy/C(or Gy/rdg)