

# Development of an Independent Audit Device for Remote Verification of 4D Radiotherapy

Jason Shoales, B.S.

David Followill, Ph.D.

Geoffrey Ibbott, Ph.D.



Work supported by PHS grant CA10953 and  
CA81647 awarded by the NCI

# Introduction

- Radiation therapy of thoracic and abdominal tumors is complicated by respiratory induced motion
- Larger margins are required around target volumes to account for intra-fractional movement
- Margin reduction can lead to decreased dose to normal tissue and allow for target dose escalation and improved local tumor control

# Introduction

- Margin reduction accomplished primarily by explicitly accounting for respiratory motion
  - Internal target volume (ITV)
  - Gating
- 4D techniques require a quality assurance device
  - Equipment acceptance and commissioning
  - Dosimetric verification of treatment plans
- The Radiological Physics Center (RPC) requires an independent audit device for 4D protocols

# Objectives

- Develop and build a quality assurance (QA) system that simulates respiratory induced motion
- Image the QA system along with the RPC anthropomorphic lung-thorax phantom
- Irradiate the QA system and assess dose delivery

# Development of a QA system to simulate respiratory motion

- Motion extents derived from literature
- Typical sinusoidal human respiratory pattern was modeled.
- Assume correlated motion between external surfaces and internal tumor motion.
- Use a two plate design for simplicity and rigidity.


# Results (Development)

- Phantom QA System Construction
  - Custom software developed for creation and storage of respiratory patterns
  - Simple, stand-alone operation with robust design suitable for travel
  - Modify original design to accommodate CT field of view size limitations




# Results (Development)

- Phantom QA System Construction



# Image QA System along with RPC Thorax-lung Phantom

- Anthropomorphic phantom containing simulated heart, spine, lung and tumor.
- Acquire static, breath-hold, free-breathing, gated and 4D CT datasets.
- Assess imaging reproducibility by measuring internal structure volumes
- Clinically constrained treatment plans






# Results (Imaging)

- Imaging Reproducibility
  - Qualitative agreement with known tumor volumes for static, BH, 4D, free-breathing and gated CT acquisitions
  - Reproducibility decreased with motion while gating demonstrated the largest range of volumes
  - Central axial diameter of the tumor shows greatest tumor shape distortion with gated datasets

# Results (Imaging)

## Tumor Volume Distribution



# Results (Treatment Planning)

- Techniques
  - Static
  - Free-breathing ITV
  - Gated
- Treatment Plans
  - Clinically constrained prescriptions
  - Three fields, 6 MV
  - 20 Gy prescription for at least 95% PTV coverage
  - Pinnacle<sup>3</sup> (v. 6.2b) used



# Dosimetric Evaluation of the Phantom QA System

- Dosimeters
  - TLD's
    - absolute dose in tumor (superior and inferior), heart, and cord.
  - Radiochromic Film
    - 2D dose distributions (axial, sagittal, coronal).
    - Isodose comparisons, binary agreement maps and profile analysis.
  - Dosimetry reproducibility evaluation
    - Static technique used as baseline, 3 trials each for free-breathing ITV and gated.
- Criteria
  - $\pm 5\%$  or 3mm distance to agreement (TG-53)

# Results (Dosimetry)

- Dosimetry Assessment
  - High level of agreement with calculated values for all TLD measurements (within 1% for all techniques at the tumor location)
  - Excellent agreement shown at the  $\pm 5\%/3\text{mm}$  level for axial, coronal and sagittal films for static and free-breathing ITV techniques.
  - Expected regions of disagreement in the beam penumbra for static and free-breathing ITV techniques.
  - Poor agreement outside the tumor volume for the gated technique


# Results (Dosimetry)

## TLD Results Summary (Dose values in cGy)

Technique	TLD	Pinnacle <sup>3</sup>	TLD/Pinnacle <sup>3</sup>
Static	2125.7	2102.4	1.011
Free-Breathing ITV	2142.7	2118.1	1.012
Gated	2114.7	2119.6	0.998

# Results (Static Plan)



Average axial AP profile from normalized film and Pinnacle<sup>3</sup> calculated profile comparison for the static plan



# Results (Static Plan)

Isodose comparison of film (dashed) vs. calculated (solid)


Binary agreement map results from 5%/3mm – 7%/7mm





# Results (FB ITV Plan)



Average axial AP profile from normalized film and Pinnacle<sup>3</sup> calculated profile comparison for the free-breathing ITV plan



# Results (FB ITV Plan)


Isodose comparison of film (dashed) vs. calculated (solid)

Binary agreement map results from 5%/3mm – 7%/7mm



# Results (Gated Plan)



Average axial AP profile from normalized film and Pinnacle<sup>3</sup> calculated profile comparison for the gated plan



# Results (Gated Plan)


Isodose comparison of film (dashed) vs. calculated (solid)

Binary agreement map results from 5%/3mm – 7%/7mm



# Results (Gated Plan)

Results for 3 film planes





Axial

Coronal

Sagittal


# Results (Gated Plan)

- Gated irradiation revealed tumor centroid displacement in the anterior and superior directions due to respiratory motion CT artifact.



# Results (Gated Plan)

Improved results for gated data with Pinnacle<sup>3</sup> data shifted.



Axial


Coronal

Sagittal

# Results

## (Binary Comparison Summary)

Overall Binary Agreement Map Histogram Analysis





# Conclusions

- The phantom QA system was simple to use and produced reliable operation throughout several trials of use.
- The phantom QA system works well in a typical clinical environment and produced good results for CT imaging reproducibility, treatment planning and treatment delivery for the static and FB ITV techniques.
- The phantom QA system indicates that the gated technique should be used with caution during the simulation and treatment planning stages to avoid a potentially serious error in radiation isocenter placement.

# Acknowledgements

- Peter Balter, Ph.D.  
George Starkschall, Ph.D.  
Naresh Tolani, M.S.  
Mark Bushman
  - M.D. Anderson Cancer Center, Houston, TX.
- Bryan Stewart
  - Radiological Physics Center, Houston, TX.

